

DOMINISH

Oxygen Bleach Powder

→ DESCRIPTION

Multipurpose powder containing an oxygen bleach for destaining, and detergents and alkalis for cleaning.

→ BENEFITS

- Colour-safe oxygen bleach
- Much safer to use than chlorine bleach
- Contains a blend of three enzymes for removal of protein, starch and mannan gums
- Detergents, enzymes and oxygen bleach work together to remove a wide range of stains
- Highly concentrated for maximum economy
- Contains 2.3% available oxygen

→ APPLICATION

Dominish is suitable for use in the laundry as a soaker and in wash booster, as well as a kitchen destainer.



→ TECHNICAL INFO

Fragrance	Colour	Form	pH Level
No added fragrance	White	Powder	10.8

Incompatible surfaces

Do not use on wool, silk or other delicate fabrics

→ AVAILABLE IN

Product	Code
8kg Pail	C12146



→ DIRECTIONS FOR USE

Dominish works best in hot water.

LAUNDRY BOOSTER

Use in addition to regular laundry powder or liquid to assist with stain removal. Increase dosage if soiling is heavy.

TOP LOADER:

WASHING MACHINE CAPACITY				
Machine Size	4kg	6kg	8kg	10kg
Grams per wash	20g (½ scoop)	40g (1 scoop)	60g (1½ scoops)	80g (2 scoops)

FRONT LOADER:

WASHING MACHINE CAPACITY				
Machine Size	4kg	6kg	8kg	10kg
Grams per wash	10g (½ scoop)	20g (½ scoop)	30g (¾ scoop)	40g (1 scoop)

LAUNDRY SOAKER

- Use at a ratio of 40g (1 scoop)/20L of hot water
- Do not use on wool, silk or other delicate fabrics

DESTAINING CROCKERY AND OTHER KITCHEN UTENSILS

- Use at a ratio of 40g (1 scoop)/10L of hot water
- Soak items for one hour
- Rinse with fresh water

Dominant
dominant.com.au

1300 789 852 enquiry@dominant.com.au

F
ALKALINE

A Safety Data Sheet for this product can be found at dominant.com.au or use the QR code



SAFETY INFORMATION

RECOMMENDED PPE

24 Hour Medical Emergency Line 13 11 26