

SaniPlus Wipes

Hospital Grade Disinfectant Wipes

→ DESCRIPTION

SaniPlus is a water-based Hospital Grade Disinfectant wipe suitable for hard surface cleaning and disinfecting within hospitals and other areas.

→ BENEFITS

- TGA listed hospital wipe AUST L 474205
- 2-in-1 cleaning and disinfection
- Proven to kill bacteria including *S.aureus*, *P.aeruginosa* and *E.hirae*
- Kills COVID-19 and Vaccinia virus
- Fragrance free
- Alcohol free
- Easy to use and convenient resealable and flexible packaging
- Non-corrosive to metals
- Durable thick spun lace fabric

→ APPLICATION

SaniPlus Hospital Grade Disinfectant cleaner is designed for use in medical settings such as hospitals, dental clinics and aged care facilities.

→ DIRECTIONS FOR USE

- Remove a wipe from SaniPlus packaging ensuring the resealable lid is closed after removal.
- Wipe dirty surface in an S-shaped action. For disinfection, pre-clean hard surface, then apply wipe onto surface.
- Allow surface to remain wet for 5 minutes and allow to air dry.
- Discard wipe after use (Do not flush in toilet).
- Use a new wipe for a new surface.

TECHNICAL INFORMATION

- SaniPlus wipes contain a unique formulation with synergistic blends of ingredients to create outstanding cleaning and disinfection performance.
- Subjected to rigorous testing and assessment by the TGA which includes EN13727, EN14476 (Vaccinia Virus) under dirty conditions.
- Do not mix with detergents as disinfection will be deactivated.
- Not suitable for use on medical devices, therapeutic goods or skin cleansing.



→ TECHNICAL INFO

Fragrance	Colour	Form	pH Level
No added fragrance	White	Solid	5.0-7.0

→ AVAILABLE IN

Product	Code
12 x 60 wipes (carton)	C29854



Dominant
dominant.com.au

1300 789 852  enquiry@dominant.com.au



A Safety Data Sheet for this product can be found at dominant.com.au or use the QR code



SAFETY INFORMATION

RECOMMENDED PPE

24 Hour Medical Emergency Line 13 11 26

Issue Date: 30/01/2025 Version: 01