VERO

Antibacterial Dishwashing Liquid

→ DESCRIPTION

An all-in-one antibacterial dishwashing detergent.

→ BENEFITS

- → Cleans and sanitises dishes in one step
- → Extensively tested with proven anti-microbial results against a range of food poisoning bacteria
- → Lemon fragrance
- → Foam level indicates sanitising efficacy levels
- → Dispensed through a Dominant dispensing system for accurate, consistent dosing
- → Featuring Dominant's proprietary Protekta® Cap, a no-touch system ensuring staff are never exposed to the product



→ APPLICATION

For use as a part of the Dominant Smart Housekeeping system

\rightarrow

DIRECTIONS FOR USE

- . Wash dishes in solution
- 2. Rinse under clean water
- 3. Allow to air dry for optimum hygiene

MANUAL DISHWASHING:

Use at a dilution of 1:500 or 2mL/L.

GENERAL CLEANING:

Use at a dilution rate of 2mL/L.

If using through a dispensing bottle, do not dilute with water before adding to the bottle as this will dilute the preservative.

Do not add any soaps, detergents or other chemicals to Vero as deactivation may occur. Never add any manual dishwashing detergents as they will completely deactivate the product.

Vero is designed for cleaning and sanitising hard surfaces, such as plates, cutlery and general dishes. Do NOT use as a sanitiser on fresh fruit and vegetables. For this application, use Dominant Food Sanitiser.

Vero has been independently tested under dirty conditions at 1:500 using the principles outlined in BS EN 1276:1997 against the following organisms:

- → Salmonella choleraesuis
- → Listeria monocytogenes

Escherichia coli Vero passed the test with >log 5 reduction after 2 minutes at 40°C. This is a kill rate > 99.999%. Testing carried out by NATA registered analytical laboratory Merieux Nutrisciences.



TECHNICAL INFO

Fragrance	Colour	Form	pH Level
Lemon	Clear Yellow	Liquid	6.5



	Product	Code
	4 x 2L Protekta® Valve Bottle	C16764
	2 x 5L Protekta® Valve Bottle	C17121
	15L Drum	C17398
	Printed 750mL Dispensing Bottle	C10433
	Flip Top Cap	C19512
	Empty 1L Labelled Dispensing Bottle	C27102
	8mL Pump to suit 5L & 1L Bottles	C13525

















SAFETY INFORMATION

RECOMMENDED PPE



